



HLT TABLETOP PADS

- Optimal bony prominence protection
- Radiolucent for minimal attenuation of fluoroscopic imaging
- Ultrasonic seam welded for fluid resistant properties
- Eliminates need for protective overlays



ARMBOARD PADS

- Soft, comfortable support of the patient's arms
- Designed with pressure management technology
- Varying pad heights to match tabletop pad thicknesses with or without X-ray tops
- Loop fastener strip runs the length of the pad



ORTHOPEDIC TABLETOP PADS

- Comfortable with significant pressure management properties
- Radiolucent for minimal attenuation of fluoroscopic imaging
- Ultrasonic seam welded for fluid resistant properties
- Hypoallergenic material