

## Radiation Protection Aprons



### REVOLUTION

Built from the ground up, Revolution Aprons offer the user the best fit of any apron on the market. The Revolution Vest and Skirt blends the highest protection with the best comfort to create a remarkable x-ray protection apron, an apron that changes what we know about x-ray protection garments.

**L103**  
Revolution Lumbar  
Vest & Skirt



**103**  
Revolution  
Vest & Skirt



**203**  
Revolution  
Full Wrap Black Belt



**303**  
Revolution  
Velcro® Front



**403**  
Revolution  
Elastic Back Front



**503**  
Revolution  
Black Belt Front

## Radiation Protection Aprons



### CLASSIC

Infab Classic Aprons have been an industry standard for over 30 years. With over 100+ in stock colors to choose from, Classic Aprons give customers the ability to choose the exact look they want. Different colors inside and out, different trim and Velcro colors and the ability to create a custom fabric, Classic Aprons are the most customizable aprons available.

VAS  
Vest & Skirt



AWC  
Weight Reliever



LSW  
Black Belt Wrap



LSF  
Black Belt Front



SDO  
Surgery Drop Off

## Lead Glasses

### Nike Frames

Infab is pleased to introduce Nike brand glasses to our Revolution line of radiation protection products. Nike is known worldwide for their high-quality products and now you can get Nike frames for your radiation protection eyewear. Available in 6 styles, every pair of Nike radiation protection glasses are comfortable, stylish and lightweight.



BANDITS  
Nike



BRAZEN  
Nike



CRUISER  
Nike



NI7085  
Nike

### Microlite Maxx 30

First protective eyewear to receive the American Ophthalmology Associations seal of acceptance. Meets ANSI Z87.1+ standards along with ASTM F803 impact resistance standards. Unisex shape and temple bars and nose bridges that have soft inside rubberized padding for extra comfort. Every pair comes complete with their own semi-rigid case and adjustable nylon





## Scatter Armor Shields



### SCATTER ARMOR

#### Scatter Protection


Scatter Armor Shields provide scatter radiation protection for radial and many other types of procedures. KIARMOR Bi-Layer is the most protective core material against scatter radiation.

#### Absorbent

Scatter Armor Shields utilize the highest absorbency in the industry, critical in heavy fluid situations.

#### How it works

Scatter Armor Shields are placed on top of the patient drape to block scatter radiation at the source, creating a scatter-free zone, reducing harmful radiation exposure to the doctor or clinician.



#### Scatter Armor Multipurpose Shield

Shunt, biliary drainage, spinal, stents, and diagnostics



#### Scatter Armor Multipurpose Shield with Fenestration

Needle biopsy or puncture procedures



#### Scatter Armor EP Shield

Subclavian access: pacemaker implants, ICS procedures, and biventricular pacing.




#### Scatter Armor Angiography Shield

Angiography, stents, coronary catheterization



#### Scatter Armor Biliary Split Shield

Biopsy procedures



#### Scatter Armor Dialysis Drape

Dialysis access and declotting procedures.



The only lead free core material which passes DIN (6875-1), ASTM (F3094-14) and IEC (61331-1-2014) standards from 50 to 150 kVp.

Disposable Thinking Cap  
DC25



## Lead Gloves

### REVOLUTION

#### Latex and Powder Free

Every pair of Revolution Radiation Protection Gloves is completely latex and powder free.

#### Tactile Grip

Revolution Radiation Protection Gloves come with tactile grip finger tips. The tactile grip helps with keeping a strong grasp on instruments during a long procedure.

#### Three Protection Levels

There are three protection levels for Infab's Revolution Radiation Protection gloves: Light, Medium, and Heavy. Choose what best fits your needs.



Type	Direct Beam Energy Level		
	60 kVp	80 kVp	100 kVp
LIGHT	40%	35%	30%
MEDIUM	68%	61%	54%
HEAVY	75%	69%	63%

Values can vary by + 5%. 2.5 mm Al - Total Filtering



**REV-TC**  
Revolution Thyroid Collar



**SVT**  
Visor Style Thyroid Collar

## Thyroid Collars

Infab thyroid collars offer a better fit and a larger area of protection. The Revolution fabric makes for a much more comfortable product that won't irritate the neck even during a long procedure. Different in look, the Revolution Thyroid collar is the only collar that gives you complete protection and total comfort.



**TCV**  
Thyroid Collar with Binding



**NST**  
Thyroid Collar without Binding

# SM Synergy Medical Corporation



REVOLUTION  
Vest & Skirt Rack



Budget Mobile  
Apron Rack  
with 10 Arms



10 Hanger  
Mobile Apron Valet



Mini Mobile  
Apron Rack



Wall Mounted Rack

## Lead Apron Racks

Infab Apron Racks are well-designed and thoughtfully engineered to provide the most efficient functionality. Each Infab Mobile Apron Rack offers locking brakes, high quality hospital-grade casters with brakes on two wheels for ease of movement and high grade, powder coated steel construction for years of trouble free service.

## Lead Barriers

Infab Barriers incorporate years of advancements and features to produce industry leading quality, innovation and protection. All barriers are available with either lead glass or acrylic. There are currently 10 different styles to choose from, making it easy to find a barrier that matches your exact specifications. If one of our 10 styles does not meet your specifications, we will work with you to create a custom barrier that can fit any desired size. With over 30 years of experience we know we can help you build exactly what you need.



Physician Protection Mobile Screen



Standard Mobile Barrier



**REVOLUTION**  
Premium Mobile Curtain  
REV-PMC





## Casette Covers

Infab cassette covers have a wide range of benefits, whether it's protecting from the spread of bacteria, or to simply keep your cassettes clean and free of debris. Preventing the spread of bacteria is essential in any medical facility; therefore the need for cassette covers is essential. Infab's cassette covers are specifically designed to fit and protect your cassettes from blood and any other potentially infectious fluids.

### Cassette Covers

Box of 100/500



### Zip Top Cassette Covers

Box of 100/500



Zip Top



Roll and Rack

## Accessories



CLX28  
Apron Cleaner



Revolution Elastic Back Thinking Cap



PGS  
Gonad Shields



REV-BAG  
Revolution Apron Bag



PTD  
Patient Transfer Device