

PORIFEROUS Polyethylene Implants Featuring SU-POR[®] Biomaterial.

PATENT PENDING

Poriferous is pleased to offer the following patent pending device



Su-Por AIRO[™] Implant

Designed with Michael A. Burnstine, M.D.

“The Su-Por AIRO Implant is a composite microchannel-suture brow suspension implant comprised of a porous sheet with a microchannel inlay. The implant acts to lift and fill the eyebrow. Where other suture and barbed suture techniques have failed because of cheese-wiring, this minimally invasive brow suspension technique succeeds because of elevation of the ROOF fat with porous scaffold to allow tissue ingrowth to help stabilize, elevate, and fix the eyebrow. The implant is anchored and stabilized to the deep temporalis fascia at the hairline for temporal eyebrow elevation.”

-Michael A. Burnstine, M.D.

Cat #	Description	Dimensions
4442	Airo -Small Qty 1	10mm x 5.0mm x 0.45mm
4443	Airo -Large Qty 1	21mm x 9.5mm x 0.45mm
4444	Airo -Small Qty 2	10mm x 5.0mm x 0.45mm
4445	Airo -Large Qty 2	21mm x 9.5mm x 0.45mm

New
Patent Pending



Scan for video
demonstration

All content herein is protected by copyright, trademarks, and other intellectual property rights owned by or licensed to Poriferous, LLC. or its affiliates unless otherwise indicated, and must not be redistributed, duplicated, or disclosed, in whole or in part, without the express written consent of Poriferous, LLC. This material is intended for health care professionals. Distribution to any other recipient is prohibited. Poriferous, LLC. does not practice medicine. The treating surgeon is responsible for determining the appropriate treatment, technique(s), and product(s) for each individual patient. Check for country product clearances and reference product specific instructions for use. For Poriferous's privacy statement, complete product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the package insert or visit our website at www.poriferous.com or call 770.683-3855. Form No. QF-171106-1 - Rev A