



Slider Style Patient Transfer Boards

Constructed of HDPE, 3/16" thick. Non-anti-static. Generous handle for easy transfers. Easily transfer patients from bed to gurney or OR Table. Smooth easy to clean surface. Available in two standard sizes.

121	22"W x 33"L x 3/16"H
122	22"W x 72"L x 3/16"H



Patient Transfer Roller Boards

Roller style patient transfer boards easily transfer patients from bed to gurney to OR table. Constructed of aluminum for easy to clean and minimal weight reasons. Standard style patient transfer board, 15" width and 3 lengths to meet your needs.

123	25"W Patient Transfer Roller Board
125	30"W Patient Transfer Roller Board
130	67"W Patient Transfer Roller Board



Wall Rack for Patient Transfer Roller Boards

Constructed of stainless steel material for ease of cleaning and prevention of infection spread. Fits all types of roller boards. Easy to hand and keeps boards clean and out of the way, yet in the spot you need them. Each rack holds one board. Mounting hardware not included due to varying substrates.

132	Wall Rack for Patient Transfer Roller Boards
-----	--



Wall Rack for Slider Style Patient Transfer Boards

Constructed of stainless steel material for ease of cleaning and prevention of infection spread. Hook holds small or large size patient transfer slider boards. Easy to hand and keeps boards clean and out of the way, yet in the spot you need them. Each rack holds one board. Mounting hardware not included due to varying substrates.

120

Wall Rack for
Slider Style Patient
Transfer Boards



Replacement Covers for Patient Transfer Roller Boards

Replacement covers for patient transfer roller boards. Fits corresponding length of roller board. Not guaranteed to fit other brands of roller boards.

124

25"W Replacement
Cover for Patient
Transfer Roller Board

126

30"W Replacement
Cover for Patient
Transfer Roller Board

131

67"W Replacement
Cover for Patient
Transfer Roller Board

